EXCELLNHP

endoMED

ELETTROBISTURI PER CHIRURGIA MONOPOLARE E BIPOLARE
ELECTROSURGICAL UNIT FOR MONOPOLAR AND BIPOLAR SURGERY
ELECTROBISTURÍ PARA CIRUGÍA MONOPOLAR Y BIPOLAR

ISO 9001
9120.ALSA

IT - 1231

ISO 13485
9124.ALS2

CE 0051

alsa apparecchi medicali

alsa apparecchi medicali

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
ADC System - "Constant Power" ADC System
 APC System - "Constant Voltage" APC System
 SIC System - "Starting Impulse Control" System
 SBC System - "Saline Bipolar Control" System

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Simultaneous activation by two surgeons
Use in monopolar endoscopy under liquid
Use in operative flexible endoscopy
Bipolar coagulation and sealing of big vessels
Bipolar cut and bipolar endoscopic cut in saline solution

The EXCELL NHP Endomed

is an electrosurgical unit complete, technologically extremely advanced and suitable for every kind of monopolar and bipolar use in all the surgical fields.

It is even provided with a special current for the bipolar cut under liquid, and it is self controlled by a specific and innovative operative software, which makes it particularly indicated for all the new endoscopic procedures in urology and gynaecology in saline solution.

10 MONOPOLAR CURRENTS

5 FOR PURE CUT AND COAGULATING CUT

- **Pure** - sinusoidal not modulated current for pure cut without coagulating effect
- **Blend 1** - sinusoidal modulated current for coagulating cut
These two currents, thanks to their specific control software, are suitable for all kinds of use, including the endoscopic procedures under liquid (TURP, TURV)
- **Blend 2** - modulated current for the cut with a strong coagulating effect Spray type for surgery or laparoscopy
- **Auto Pure** - sinusoidal not modulated current for cut without coagulating effect
- **Auto Blend** - sinusoidal modulated current for coagulating cut

1 FOR CUT IN FLEXIBLE ENDOSCOPY

- **Auto Endo** - current with alternating phases of cut and coagulation for the flexible endoscopy

4 FOR COAGULATION

- **Fulg Forced** - modulated high voltage current with strong coagulating superficial and deep effect
- **PinPoint-Contact** - modulated medium voltage current with medium coagulating superficial and deep effect
- **Soft** - modulated low voltage current with delicate coagulating effect without superficial carbonization
- **Spray** - modulated very high voltage current for a very strong superficial coagulation with a low penetration into the tissues, even without any contact of the active electrode

100 MEMORY PROGRAMS

The Excell NHP Endomed allows memorizing and then selecting, by pushing a simple button, 100 complete function programs which can be memorized by the operators. The memories include pre-set programs for the general or the laparoscopic surgery, the monopolar endoscopy under liquid, the operative flexible endoscopy, as well as the bipolar mode or the bipolar endoscopy in saline solution.

The unit memorizes the last adjustments used when switched off.

NEUTRAL ELECTRODE CONTROL "NPCC SYSTEM"

The circuit can be used with all kinds of neutral electrode, "Non Split - single section" or "Split - dual section" type.

Using "Non Split" electrodes, it checks the cable and its connections.

Using "Split" electrodes, it verifies also the contact between the electrode and the patient tissues.

It stops the power delivery, with acoustic and visual alarm, if:

- The connecting cable gets broken or disconnected.
- The contact gets insufficient (Using "Split" electrodes only).

HF LEAKAGE CONTROL

A specific circuit continuously checks the high frequency leakage currents to earth. In case of potentially dangerous situations, like the contact between the patient and the metallic parts of the OT-table, it automatically decreases the output powers to avoid the risks involved (i.e burnings of the patient), and it informs the operators by an alarm signal.

TWIN FOOT-SWITCH

One twin foot-switch, selectable for monopolar or bipolar mode, protected against the involuntary activation, waterproof (IPX7) and with medical safety voltage supply (12VDC).

2 MONOPOLAR OUTPUTS

The simultaneous activation of the two outputs by two surgeons is always possible.

The outputs can be used both by hand-switch handles, with International 3-pin connector, or one by hand-switch handle and one, by the twin foot-switch, using a normal monopolar handle or an instrument cable.

To use all the monopolar instrument cables, with single connector different from the standard one, the adaptor RD/5 is available.

MAXIMUM CONTROL MAXIMUM SAFETY

The auto-check system "Master/Sleeve System with two Microprocessors at 32 bit with interconnection in feed-back" grants maximum functionality performing:

- **Main Auto-Check** - The complete hardware/software auto-check at switching on, then repeated every 40 minutes;
- **Standard Auto-Check** - The functioning control during the use, even including specific controls about the outputs or the activation ways, acceptable according to the selected modes.

The self-check system stops the power delivery, informing the users by giving specific Error Codes, if it detects problems about:

- **General Error Control** - The general functioning, the accessories connection, the usage or the activation;
- **Output Power Control** - The powers delivery;
- The system memorizes the last **32 Error codes** to allow further checks.

COMPLETE "OUTPUT POWER CONTROL"

All the currents are controlled by microprocessor through four different software systems:

- **ADC System "Constant power"** which, according to a continuous feed-back between the unit and the patient tissues, automatically optimizes the power delivery in order to intervene in a very efficacious way – with the lowest possible power – on many tissues or under liquid;
- **APC System "Constant voltage"** which, according to a continuous feed-back between the unit and the patient tissues, automatically adjusts the output powers, keeping constant the currents electric arc, in order to intervene on the tissues, when necessary, in a delicate way with the minimum superficial necrosis;
- **SIC System "Starting impulse control"** which automatically checks the currents starting impulse in order to grant their maximum efficacy using the lowest possible power;
- **SBC System "Saline Bipolar Control"** which, according to a continuous feed-back between the unit and the operating field, automatically optimizes the efficacy of the bipolar cut under liquid in saline solution.

MANUFACTURING TECHNOLOGY ON-SITE CALIBRATION SOFTWARE UPGRADE

The manufacturing technology grants the highest reliability. The HF currents generator is "**Mosfet powered damped oscillator with switching supply**" type for use without time limit. To help technical maintenance, the Excell NHP Endomed is made of interchangeable boards and it memorizes the last it **32 Error Codes**; it can be on-site calibrated and moreover the software updating can be made by a serial port connected to a PC.

SMOKES EVACUATION AND FILTERING

The Excell NHP Endomed can be completed by automatic suction units that grant evacuation and safe filtering of the surgical smokes, by reducing the risks of contamination and by solving the problem of the visibility during the laparoscopic procedures.

5 BIPOLAR CURRENTS

1 FOR CUT

1 FOR ENDOSCOPIC CUT IN SALINE SOLUTION

2 FOR COAGULATION

1 FOR SEALING OF BIG VESSELS

- **Pure** - pulsed current for normal cut
- **Saline** - pulsed current for endoscopic cut in saline solution
- **Micro** - pulsed current for coagulation Soft / Micro Precise type and for coagulation in saline solution
- **Macro** - pulsed current for coagulation Standard / Forced type in laparoscopy
- **Seal** - pulsed current with automatic stop for sealing of big vessels until diameter of 7 mm

OUTPUT CIRCUIT

The output circuit is "Floating - insulated at high and low frequencies type" and protected against the use of the defibrillator. It is the best one to grant the safest use. The IEC 601-1 safety classification of the unit is CF, the highest and the only one allowed for direct application on the heart.

1 BIPOLAR OUTPUT

The output can be activated by the twin foot-switch or, only selecting the Seal current for sealing of big vessels, by the foot-switch with automatic stop. To connect all the instrument cables, with connector different from the ALSA standard, specific adaptors are available.

MONOPOLAR AND BIPOLAR CURRENTS

MONOPOLAR CUT

Currents	Currents Features					Currents self regulation			
Pure	350 W	350 Ω	3450 Vpp	CF: 1.6	M: no DT: no	ADC System		SIC System	
Blend I	300 W	350 Ω	3600 Vpp	CF: 2.3	M: 29 kHz DT: 65%	ADC System		SIC System	
Blend II	140 W	600 Ω	7600 Vpp	CF: 8.1	M: 19 kHz DT: 9%	ADC System		SIC System	
Auto Pure	350 W	350 Ω	1350 Vpp	CF: 1.6	M: no DT: no		APC System	SIC System	
Auto Blend	300 W	350 Ω	1930 Vpp	CF: 2.3	M: 29 kHz DT: 65%		APC System	SIC System	
Auto Endo	220 W	350 Ω	1710 Vpp	CF: 2.2	50% Pure 50% Blend I		APC System	SIC System	

MONOPOLAR COAGULATION

Currents	Currents Features					Currents self regulation			
Fulg Forced	150 W	350 Ω	4700 Vpp	CF: 4.5	M: 78 kHz DT: 3,5%	ADC System			
PinPoint Contact	250 W	250 Ω	3460 Vpp	CF: 2.6	M: 29 kHz DT: 50%	ADC System			
Soft	280 W	250 Ω	3440 Vpp	CF: 2.5	M: 29 kHz DT: 56%	ADC System			
Spray	140 W	600 Ω	7600 Vpp	CF: 8.1	M: 19 kHz DT: 9%	ADC System			

BIPOLAR CUT

Currents	Currents Features					Currents self regulation			
Pure	180 W	350 Ω	1200 Vpp	CF: 1.5	M: no DT: no	ADC System		SIC System	
Saline	320 W	50 Ω	1200 Vpp	CF: 1.5	M: no DT: no	ADC System		SIC System	SBC System

BIPOLAR COAGULATION

Currents	Currents Features					Currents self regulation			
Micro	130 W	100 Ω	420 Vpp	CF: 1.7	M: no DT: no		APC System	SIC System	
Macro	130 W	200 Ω	1050 Vpp	CF: 1.7	M: no DT: no	ADC System		SIC System	
Seal	130 W (200 W)	100 Ω	420 Vpp	CF: 1.7	M: no DT: no		APC System	SIC System	

Currents Features

Currents self regulation electronic systems

W Output powers
Ω Nominal Load
Vpp Open circuit Peak to Peak voltage
CF Crest Factor
M Modulation
DT Duty Cycle
(W) Starting Impulse

ADC System - "Constant Power" ADC System
APC System - "Constant Voltage" APC System
SIC System - "Starting Impulse Control" System
SBC System - "Saline Bipolar Control" System

TECHNICAL FEATURES

Standards Applied	Unit is compliant with IEC 60601-1, IEC 60601-1-2 and IEC 60601-2-2 Standards.
Class and EC Approval	Class IIB – CE Approval (IMQ - CE0051): 187/MDD. 93/42EEC Certificate and subsequent amendments and additions (2007/47/EC).
Class and Type IEC 60601-1	Class I -Type CF, for use on the heart (LF leakage currents on the patient < 0.004mA).
Output Circuit IEC 60601-2-2	Floating-insulated at high and low frequencies (High Frequency leakage currents: 100mA.) Protected against the use of the defibrillator.
EMC Compatibility	Unit is compliant with IEC 60601-1-2 – Test report: IMQ nr.80SF00709/1, 80SF00709/2.
Memories and Programs	100 memorisable programs. Pre-set programs for general / laparoscopic surgery, monopolar or bipolar endoscopic use under liquid, flexible endoscopic surgery, bipolar use only.
Functioning Control	In real time by two microprocessors “Master Sleeve System at 32 bit”.
Auto Diagnosis - Faults Control - Output Error Control - HF Leakage Current Control - Error Codes Memorization -	Complete hardware/software self control with: – Main Auto-Check at the switching on, then repeated every 40 minutes; – Continuous Standard Auto-Check during the whole functioning. Interruption of the functioning (in less than 100millisec) with alarm signal to the operators through specific Error Codes, when the auto-diagnosis systems find out some problems about: – General Error Control - The general functioning or the activation mistakes; – Output Error Control – The delivery of the powers. HF Leakage Control – Continuous control, by a specific circuit, of the HF leakage currents to earth with automatic decreasing of the output powers and alarm signal, just in case of dangerous situations like the contact between the patient and the metallic part of the operating table. Memorization of the 32 last Error Codes in order to facilitate further checks.
Output Powers Automatic Adjustment	Output Power Control, complete automatic regulation by microprocessors of the output powers through four different modes: – ADC System – Constant power, which auto-adjusts the powers by verifying the voltage and the current, according to a feed-back in real time (7,000 checks/sec) between the unit and the patient tissues; – APC System – Constant voltage, which auto-adjusts the powers, by keeping the voltage constant, according to a feed-back in real time (7,000 checks/sec) between the unit and the patient tissues; – SIC “Starting Impulse Control”, which auto-adjusts the starting pulse of the currents; – SBC “Saline Bipolar Control” which auto-adjusts, according to a feed-back in real time (7,000 checks/sec) between the unit and the working area, the current for the bipolar cut in saline solution.
Neutral Electrode Safety Circuit	NPCC System for use of “Non Split electrodes with one section” and “Split electrodes with double section”. Using the “Non Split” electrodes, it checks the cable and its connections (with the unit / the electrode). Using the “Split” electrodes, it checks not only the cable or its connections, but even the contact between the electrode and the patient tissues. If the cable is broken / non connected or (only using the “Split” electrodes) if the contact is dangerous, the circuit interrupts the delivery of the power and gives a luminous acoustic alarm (red light).
Monopolar / Bipolar Sockets - Activation Modes - Setting signals -	Two monopolar outputs, which can be used at the same time by two operators: – Both of them by the hand switch handles; – One by hand switches and one by the double pedal foot-switch. One bipolar output, which can be always activated independently from the monopolar ones: – For the cut or the coagulation, by the double pedal foot-switch; – Using the current for sealing of big vessels, by the pedal foot-switch with Auto Stop. Luminous signals to point out the selected and usable sockets or activation modes.
Pedal foot-switch	One double pedal foot-switch, selectable for the monopolar or bipolar functions, protected against the involuntary activation, waterproof (IPX7) and with medical safety voltage supply (12VDC).
Electronic Generator - Working Frequency - Cooling - Mains Supply - Mains Absorption at 230V - Supply Cable - Protective Earth Connection - Equipotential Conductor -	“Mosfet powered damped oscillator with switching supply” for use without time limits. Monopolar and Bipolar 440 kHz±10% By convection, without fan. Volt: 230/220-115/110 ±10% - Hz: 50-60 - Phases: alternated single phase. Max power: 3.6 A = 828 VA - At switching on: 3.65 A = 840 VA - Stand-by: 0.4 A = 92 VA. Detachable type with SHUKO plug. Included in the supply cable. Resistance < 0,1Ω. International standard plug DIN42801.
Currents Setting	By push-buttons with settings pointed out by luminous signals.
Powers Adjustment	Micro and Macro by push-buttons with control by display.
Colour Codes and Signals	Colour codes and luminous or acoustic signals of activation and alarm compliant with IEC 60601-2-2.
Enclosure and Control Panel	Metallic enclosure protected against inlet of liquids according to IEC 60601-2-2. Flat panel.
Dimensions and Weight	L x D x H: cm 38x38x16 - Kg 15.
Technical Assistance - Technological Updating - Calibration -	Memorization of the latest 32 Error Codes in order to facilitate the check and the technical assistance. Software updating through a serial port connected to a PC. On - Site calibration by specific software.
Use and Keeping	Temp 10°-40°C / Humid. 35-75% / Press. 700-1060 hPa (Standards according IEC 60601-1).

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UNIT AND STANDARD SETS OF ACCESSORIES

B960E - EXCELL NHP ENDOMED, without accessories

B610/A - SET OF STANDARD ACCESSORIES composed by:

- DS/E** Twin foot-switch (x1)
- NP/A** Stainless steel neutral electrode for adults - cable 2.5mt (x1)
- FGE** Fixing belt for neutral electrodes with 2 buttons (x1)
- MPE/E** Sterilizable electrodes-holder handle, cable 3.5mt (x2)

SEL/E - SERIES OF ELECTRODES, composed by:

- E1** Straight knife electrode (x2)
- E5** Thick needle electrode (x2)
- E7** Thin needle electrode (x2)
- E12** Ball electrode Ø 2.5mm (x1)
- E14** Ball electrode Ø 4mm (x1)

B610/B - SET OF STANDARD ACCESSORIES same as **B610/A**, but with flexible neutral electrode made of conductive rubber for adults **NP/GA**

B610/P - SET OF STANDARD ACCESSORIES same as **B610/B**, but with paediatric neutral electrode **NP/GP**.

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NEUTRAL ELECTRODES

- NP/A Stainless steel (13x25cm), cable 2.5mt
- NP/GA Flexible conductive rubber (12x25cm), cable 2.5mt
- NP/GP Flexible conductive rubber (8x18cm), cable 2.5mt
- FGE Fixing belt for neutral electrodes with 2 buttons (*)

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DISPOSABLE, ADHESIVE NEUTRAL ELECTRODES

- CMS/E Re-usable connecting cable 2.5mt long
- EIP/SA Adhesive, Split double section electrode, adults type (25pcs in each pack.) (*)
- EIP/SP Adhesive, Split double section electrode, paediatric type (25pcs in each pack.) (*)

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Electrodes for microsurgery, sterilizable

MID	Reducer for needles (for all the electrodes-holders)
AID	Insulated needle Ø 0.10mm
SAD	Series of 12 needles Ø 0.10mm
SAD/1	Ditto, Ø 0.15mm
SAD/2	Ditto, Ø 0.20mm
SAD/3	Ditto, Ø 0.40mm

- **Pinze, strumenti e cavi monopolari per chirurgia, endoscopia e laparoscopia.** *Vedi specifico catalogo.*
- **Monopolar forceps, instruments and cables for surgery, endoscopy and laparoscopy.** *See the specific catalogue.*
- **Pinzas, instrumentos y cables monopolares para cirugía, endoscopia y laparoscopia.** *Ver catálogo específico.*

- **Resettoscopio bipolare con accessori.** *Vedi specifico catalogo.*
- **Bipolar resectoscope with accessories.** *See the specific catalogue.*
- **Resector bipolar con accesorios.** *Ver catálogo específico.*

- **Forbici bipolari per chirurgia, Clamps per chiusura grossi vasi in chirurgia e laparoscopia.** *Vedi specifico catalogo.*
- **Bipolar forceps for surgery, Clamps for the sealing of thick vessels in surgery and laparoscopy.** *See the specific catalogue.*
- **Tijeras bipolares para cirugía, Clamps para el cierre de gruesos vasos en cirugía y laparoscopia.** *Ver catálogo específico.*

- **Pinze, strumenti e cavi bipolari per chirurgia e laparoscopia.** *Vedi specifico catalogo.*
- **Bipolar forceps, instruments and cables for surgery and laparoscopy.** *See the specific catalogue.*
- **Pinzas, instrumentos y cables bipolares para cirugía y laparoscopia.** *Ver catálogo específico.*

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