ESU-GA/MPC

Argon gas unit

Quality System Certification

ISO 9001/EN 46001

CE 0051

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ESU-GA/MPC is an "Add-on module" for Argon-enhanced electro-surgical procedures of cut and coagulation for open, laparoscopic and endoscopic surgery. Argon is an inert gas that does not combine with other chemical elements and, used thanks to specifically designed accessories in combination with the HF currents generated by an electro-surgical unit, allows to have 2 different, otherwise unobtainable, results:

1) The first one is a better surgical end-result because of the displacement of oxygen from the operation area eliminating the odours, smokes and granting an excellent view.

2) The second one, more interesting, is a kind of coagulation completely different from the normal ones obtained by the sole use of the electro-surgical unit. As matter of fact, when Argon gas is employed, its flow is ionized by the HF current that flows freely in such channel reaching the target tissue and generating an absolute constant superficial coagulation with a maximum depth of 3mm. Furthermore, since this flow basically moves from high-impedance tissue towards low-impedance tissue, the haemostatic effect is maximum in the non-coagulated tissues and absent in the already coagulated ones in consequence of reducing the carbonization at the maximum.

For the above described reasons, this coagulation results extremely effective and useful not only in traditional open sky surgery, i.e. for the parenchymal tissues, but also in the laparoscopic and endoscopic ones for which are designed specific electrodes. For the endoscopic surgery it is available even a very useful twin double pedal switch for activation of, both electro-surgical unit and ESU-GA/MPC.

The module it is easily connectable to the electro-surgical units by simple interconnection system, can be activated either by hand-switch handle or pedal switch and the system ESU-GA/MPC/electro-surgical unit can be operated, even simultaneously by two surgeons, both for the Argon-enhanced coagulation and the standard electro-surgical performances.

For the endoscopy has not only the facility to contemporary connect both the electrodes for polipectomy/papillotomy and the flexible probes for Argon-enhanced coagulation, but even a twin double pedal switch to activate the complete system module/electro-surgical unit.

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TECHNICAL FEATURES

- **"MPC System"** Complete control by microprocessor with specific self-test at switching on the unit and following continuous control during all the functioning with automatic stop in case of faults (**Alarm/error codes** for the users)
- **"Constant flow system and Safety gas system"** Control of the gas circuit with double safety by checking the pressure/gas input (insufficient or superior to the maximum one acceptable by the system) and by measuring of, both, the pressure of gas input (bar) and of the real flow rate at the activation (l/min).
- Interface and activation system for electro-surgical units.
- Two independent displays for gas pressure and real flow rates of gas (cut and coagulation).
- Independent manual controls for 0.2 to 12 l/min of the gas flow rates for cut or coagulation.
- Gas supply with cylinder or centralized system (max.pressure: 4 atm).
- Connector for Argon output with bacterial filter (0.2 micron) to grant its sterility.
- Activation by hand-switch handle or by double pedal switch.
- Simultaneous use of, both, Argon-enhanced coagulation and electro-surgical unit by two surgeons
- Mains/absorption/equipotential connection: 230V-50/60Hz / 20 VA / DIN connector.
- Cooling by convection, without fan.
- Hand-switch handle, double pedal-switch and active rigid electrodes for open surgery and laparoscopic use. Connecting cable, active flexible probes and twin double pedal-switch for endoscopy.
- Dimensions: cm. 35x38x12-Weight: 10 kgs
- Safety rules: Unit meets IEC 601-1 and IEC 60601-2-2 III Ed.
- Classification according to IEC 60601-1: Class I, type CF with "floating" output circuit.
- Classification and Certification according to EEC 93/42: Class IIB-CE 0051
- ALSA Quality System: Approved according to ISO 9001, VISION 2000/EN 46001

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UNIT AND ACCESSORIES

CG100 - Module ESU-GA/MPC, without accessories

STANDARD ACCESSORIES FOR OPEN AND LAPAROSCOPIC SURGERY

Using all following accessories the electro-surgical unit connected to the module needs one of its standard sets of accessories.

For open surgery with hand-switch activation of ESU-GA/MPC

CG101 - Standard set including (not included gas cylinder/pressure reducer):
AC/HANDLE - Hand-switch handle, cable mt. 3,5
AC/E25-C - Rigid electrode for coagulation, l. mm. 25
AC/E40-L - Rigid knife electrode, l. mm. 40
ESU/TG - Tube for gas input (length mt. 1 with quick - action coupling)
ESU/CP - Interface cable (ESU activation)
ESU/CA - Interface cable (HF current)

For other active electrodes for open sky surgery, please see:

General listing of accessories - AC/E electrodes

For laparoscopic surgery:

The set CG101 must be completed by AC/E 320-C electrode and ESU-CA/MAN interface cable

For open or laparoscopic surgery with pedal-switch activation of the ESU-GA/MPC:

The set CG101 must be completed by DS/E pedal.

STANDARD ACCESSORIES FOR ENDOSCOPY (GASTROENTEROLOGY)

The Electro-surgical unit used with the module requires just one neutral plate

CG102 - Standard set of accessories with pedal-switch activation of, both, ESU-GA/MPC and electro-surgical unit, including (not included gas cylinder/pressure reducer/active flexible probes):
AC/CABLE - Connecting cable for flexible probes, l. mt. 2,5
DS/T - Twin double pedal-switch, waterproof (for activation of ESU/ESU-GA/MPC)
ESU/CA-MAN - Interface cable (ESU activation and HF current)
RD5 - Adapter for ESU (for connecting HF cables of the electrodes for polipectomy, etc.)

For active flexible probes, please see:

General listing of accessories - AC/FP electrodes.

GENERAL LISTING OF ACCESSORIES

Handle, cable, electrodes - autoclavable type at 134°

AC/HANDLE - Hand-switch handle with cable (Argon gas and HF current), l. mt. 3,5
AC/E25-C - Rigid electrode for coagulation, l. mm. 25
AC/E100-C - Ditto, l. mm. 100
AC/E320-C - Ditto, l. mm. 320 (for laparoscopy)
AC/E40-A - Rigid needle electrode, l. mm. 40
AC/E115-A - Ditto, l. mm. 115
AC/E40-L - Rigid knife electrode, l. mm. 40
AC/E115-L - Ditto, l. mm. 115
AC/CABLE - Connecting cable (Argon gas and HF current) for flexible probes, l. mt. 2,5
AC/FP1 - Flexible probe for coagulation, l. mt. 1,5/diam. mm. 1,5
AC/FP2 - Ditto, l. mt. 1/diam. mm. 2,3
AC/FP3 - Ditto, l. mt. 2,2/diam. mm. 2,3
AC/FP4 - Ditto, l. mt. 2,2/diam. mm. 3,2

Trolleys, cylinder, tubes, pedal switches

H/GA - Trolley for ESU, module and 2 cylinders of 5 lt. with mains supply sockets - HxLxD: cm. 100x41x44
H10/GA - Ditto, economic type, for 2 cylinders of 5 lt. or 10 lt. without mains supply sockets - HxLxD: cm. 100x45x70
B/5 - Gas cylinder lt.5 (with gas)
RD/P - Pressure reducer for B/5 (200 bar/1 - 10 bar)
ESU/TG - Tube for gas input (length mt.1 with quick - action coupling)
DS/E - Double pedal switch, waterproof (for activation of ESU-GA/MPC)
DS/T - Twin double pedal switch, waterproof (for activation of ESU/ESU-GA/MPC)

Interface cables to the HF units

ESU/CA - Interface cable (HF current)
ESU/CA-MAN - Interface cable (HF current/ESU activation)
ESU/CP - Interface cable (ESU activation)
EEX/AC - Connector for ESU (automatic activation of setting "spray coag.")

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